REUSABLE NON-THERMAL BIOPSY FORCEPS

(ENGLISH & FRENCH IFU)

FOR COLONOSCOPES

**BX 100, BX 100N, BX 110, BX 110N, BX 140,
BX 140N, BX 910N,**

FOR GASTROSCOPES

**BX 200, BX 200N, BX 210, BX 210N, BX 240,
BX 240N, BX 300, BX 300N, BX 420, BX 520,**

FOR SIGMOIDOSCOPES

BX 250, BX 250N, BX 260, BX 260N, BX 280,

FOR DUODENOSCOPES

BX 430

FOR BRONCHOSCOPES

BX 410, BX 510.

Instructions for Use 2

Mode d'emploi 6

PRIMED INSTRUMENTS INC.
1080 Tristar Drive, Unit 14
Mississauga, Ontario, L5T 1P1
Canada
PHONE: 905-565-0565
FAX: 905-565-0566
www.primedENDO.com

EC REP

ADVENA LIMITED
Tower Business Centre
2nd Flr., Tower Street
Swatar, BKR 4013
Malta

REUSABLE BIOPSY FORCEPS

RX ONLY

CE

INTENDED USE

Biopsy forceps are used to sample lesions of the gastrointestinal tract and/or lungs, such as polyps and suspected malignancies. The physician is responsible for procedural preference and sample preparation, since these usually vary greatly. Do not use this device for any purpose other than the stated intended use.

SIGNAL SYMBOLS

WARNING	Indicates a potentially hazardous situation, which could result in death or serious injury.
CAUTION	Indicates a potentially hazardous situation, which may result in minor or moderate injury.
NOTE	Indicates additional helpful information.

POTENTIAL COMPLICATIONS

Potential complications associated with gastrointestinal endoscopy include, but are not limited to: perforation, hemorrhage, aspiration, fever, infection, allergic reaction to medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or cardiac arrest. Potential complications associated with bronchoscopy include, but are not limited to: pneumothorax, epistaxis, vocal cord injury, hypoxia, damage to teeth (from rigid bronchoscopy).

CONTRAINdications

Contraindications include those specific to primary endoscopic tissue acquisition.

CAUTIONS

Device must be sterilized prior to first use following the parameters included in this instruction. Prior to each subsequent use, the device must be cleaned and sterilized following the cleaning and sterilization instructions provided in this booklet.

Do not coil biopsy forceps into a diameter smaller than 20cm.

Do not use excessive pressure when actuating the forceps. Excessive pressure will cause damage to the instruments.

Keep forceps jaws closed during insertion into, traversal through, and removal from the endoscope channel.

Confirm compatibility with the endoscope prior to use. Using forceps with incompatible channel size will damage the endoscope.

Special Caution should be taken when handling instruments that have a needle.

When using detergent solution, always follow the detergent solution manufacturer's recommendations regarding temperature, soaking time and concentration. Always use neutral pH solution for cleaning. Subjecting biopsy forceps to highly acidic or alkaline solution for a prolonged period of time can lead to damages.

WARNINGS

Biopsy forceps should only be used by, or under the supervision of, physicians trained in endoscopic tissue biopsy.

Reusable biopsy forceps should only be used to biopsy tissue where possible hemorrhage will not present a danger to patients. Adequate plans for management of potential bleeding and appropriate airway management should be in place.

Please follow inspection instructions before using the instruments. Failure to follow instructions may lead to serious medical consequences.

Cleaning with automated reprocessing machine is NOT recommended.

Sterilization must be performed via an autoclave. High level chemical disinfection is NOT recommended.

DO NOT sterilize biopsy forceps in its original packaging. Remove biopsy forceps from original packaging and place in a suitable container/package prior to sterilization.

LIMITATIONS ON REPROCESSING

Repeated processing has minimal effect on the biopsy forceps. End of life is normally determined by wear and damage due to repeated use and/or rough handling. Follow inspection instructions closely to identify biopsy forceps that no longer can fulfill its intended use.

EQUIPMENT REQUIRED

1. Endoscope.
2. Fixative & biopsy sample container.
3. Cleaning: Ultrasonic cleaner, container, enzymatic detergent, cleaning brush, sponge, lint-free cloth.

INSPECTION

1. Remove the forceps from the package and carefully unwind the sheath from the proximal (handle) end. Remove the protective guard from the jaws.
2. Hold the distal (jaws) end of the forceps in one hand and gently operate the jaws by working the finger spool. Normally, forceps should open and close easily with two loops of approximately 20cm in diameter. **CAUTION:** Operating the forceps while the unit is still tightly coiled may lead to its damage.
3. Inspect for any functional abnormalities and any kinks or bends in the outer sheath. **CAUTION:** If any irregularities are noted, DO NOT USE the instrument. Contact PriMed Instruments Inc. or the authorized distributor for a return authorization. Do not repair instruments. Attempts to repair nonfunctional or damaged instruments may induce dangerous failure modes during use.
4. Used or reprocessed biopsy forceps that fails indicates that the device has reached its end of life. Forceps past its useful life should be treated as hazardous waste and disposed of in accordance with hospital/clinic protocol.

USE

1. First, verify that the channel size of the endoscope is adequate to accommodate the instrument with the jaws closed. The forceps jaws should be closed when introduced into the endoscope channel. Advance the instrument slowly down the channel, using small strokes, until the jaws can be visualized at the end of the endoscope. NOTE: It may be necessary to reduce the angle of the angulated portion of the endoscope to allow for ease of passage of the forceps. Do not force forceps through the endoscopes' channel.
2. To obtain a biopsy specimen, gently open the jaws and advance the forceps to the biopsy site. If needed, articulate the end of the endoscope to direct the forceps to the biopsy site. Close the jaws completely to cut cleanly through the tissues and gently pull the forceps back. NOTE: Do not apply excessive pressure when closing forceps jaws. Applying excessive force to the handle spool or thumb ring may damage the instrument.
3. Keep the jaws closed to retain the sample. Slowly withdraw the forceps. NOTE: If the biopsy sample is too large to close the jaws sufficiently, withdraw the forceps as far as the tip of the endoscope, and simultaneously withdraw the endoscope and biopsy forceps from the patient.

CLEANING

CAUTION: Immediately after use the biopsy forceps must be cleaned to remove mucus, tissue and blood from the instruments. Do not let used forceps dry before cleaning. Do not allow forceps to rest unprocessed for longer than 30 minutes. Failure to properly clean the device will compromise sterilization and may cause the device to fail. Inadequate cleaning of the forceps will void the warranty, since accumulation of dried material in the working mechanisms of the forceps will cause product deterioration.

1. Wear rubber gloves while handling contaminated biopsy forceps and working with harsh cleaning solutions.
2. Prepare a large container of warm tap water and enzymatic detergent (e.g. Enzol[®]) solution prepared according to manufacturer's instruction. The container should be large enough so that the instrument can be fully immersed, without creating coils of less than 20cm in diameter. Tighter coils may damage the device. Always loop the instrument starting at the distal (jaws) end. Place the forceps into the solution, distal end first.
3. Allow the forceps to soak for a minimum of one (1) minute. Using a soft cleaning brush or a soft toothbrush, scrub the jaws and hinges of the jaws gently but thoroughly to remove all material. Scrub the flexible sheath, if necessary, to remove any debris in the coil. If there is no debris on the shaft or handle, simply wipe it down with a sponge or gauze in the enzymatic solution.

4. Rinse the instrument in a large container with ample amount of warm tap water for minimum of one (1) minute.

5. The unit should now be placed in an ultrasonic cleaner. Fill ultrasonic tank with warm tap water and enzymatic detergent (e.g. Enzol[®]) solution prepared according to the enzymatic detergent manufacturers' instruction. Completely immerse the biopsy forceps in the ultrasonic cleaner tank. Do not exceed the 20cm minimum coils diameter. Clean the biopsy forceps in an ultrasonic cleaner with a minimum of 40 kHz frequency output for 10 minutes. Rinse thoroughly under ample amount of clean, running tap water for at least one (1) minute to remove any residual cleaning solution. Wipe dry with clean, dry, lint-free cloth.
6. Lightly lubricate the forceps distal hinges with medical grade silicone lubricant. Open and close the jaws a few times. NOTE: Proper lubrication of the forceps after each use will help maintain the performance characteristics of the device.

STERILIZATION:

WARNING: CURRENT INFECTION CONTROL GUIDELINES REQUIRE THAT BIOPSY FORCEPS AND SIMILAR ACCESSORIES THAT BREAK THE MUCOSAL BARRIER OR ENTER BODY CAVITIES MUST BE STERILIZED PRIOR TO USE.

Thorough manual cleaning prior to sterilization is essential to achieve sterility. Organic debris on the forceps and residual moisture will inhibit sterilization.

Place the biopsy forceps in a suitable container or package. Do not exceed the 20cm minimum coils diameter.

Recommended autoclave cycles:

Sterilizer	Temperature	Pressure	Exposure Time
Prevacuum	134°C	.2Mpa	4 minutes
Gravity Displacement	134°C	.310Mpa	15 minutes

STORAGE

Biopsy forceps should be stored in the sterile package at room temperature in a clean and dry area.

PRODUCT DISPOSAL

Used biopsy forceps may be a potential biohazard. Handle and dispose of in accordance with the hospital/clinic protocol and applicable regulations.

SYMBOLS MEANING

CE CE Mark. Quantity of Units in Case.

UPN Product Number. Quantity of Units in Dispenser.

This product is not made with LATEX.

This product is not made with Phthalate.

R ONLY Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

Visit us online at www.primedENDO.com for full list of symbol glossary.

NOTE: Contact PriMed Instruments if you have any questions or comments regarding the content of this document.

If two years have elapsed between the issued date for these instructions and product use, the user should contact PriMed Instruments Inc. to determine if additional/new information is available.

Enzol[®] is a registered trademark of Johnson & Johnson Inc.

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Phone: 905-565-0565
Fax: 905-565-0566
Email: info@primedENDO.com

TYPE: INSTRUCTIONS FOR USE
ID: W-7.5.3-W2C

ADDENDUM to IFU W-7.5.3-W2
**for MultiCROC® Multi-Sampling Reusable Biopsy Forceps (aka BX
MSRBF) models.**

IMPORTANT: READ AND FOLLOW ALL THE INSTRUCTIONS IN THE IFU BOOKLET W-7.5.3-W2 AND THE INSTRUCTIONS IN THIS ADDENDUM BEFORE USING THE BX MRBF.

INTENDED USE

The BX Multiple Sample Biopsy Forceps (only model numbers starting with “BX9...”) have the capability to collect multiple biopsies in a single pass of the device through the channel of the endoscope (up to 5 tissue samples per single pass). Total number of biopsies collected will depend on tissue type and user’s preference.

CAUTIONS

- Do not attempt to pass the instrument when the scope is excessively articulated. Doing so will damage the endoscope
- BX MRBF biopsy forceps jaws are designed to facilitate collection and storage of multiple biopsies. During each subsequent biopsy, tissue samples are pushed toward the back to fill-up the jaw’s elongated cups. Always monitor amount of samples stored inside the jaw’s cups and do not overload them. Overloading may result in crush artifacts and lost samples.

WARNING

- BX MRBF jaws have a smaller aperture, when opened, and a longer distal end than standard biopsy forceps. Do not apply excessive force to advance the forceps into the tissue. Applying excessive force may result in perforation of the biopsy site.
- The jaws are designed to travel about 1-2mm back during closing of the device. If necessary, compensate for this movement by gently forwarding the biopsy forceps toward the biopsy site. Make sure to apply a sufficient force to close the jaws completely. Maintain pressure to keep the jaws closed to retain the sample(s) until the next biopsy is taken.

USE

- To obtain subsequent samples, it is not necessary to withdraw the forceps from the endoscope. Up to five biopsy samples may be taken without withdrawing the forceps from the scope.
NOTE: When performing biopsy in confined spaces such as in the Esophagus, pull jaws closer to tip of endoscope to avoid potential obstruction to jaws during opening. The Jaws may not open fully unless a slight continuous pressure is applied to the handle prior to taking the sample.
NOTE: Do not drag open jaws loaded with harvested samples against the surrounding organ cavity and do not open or close the jaws rapidly. Doing so may result in a loss of tissue samples.
NOTE: If a device is to be reused on the same patient during the same procedure, be sure that all biopsy samples are removed prior to reuse. Failure to do so may result in a misdiagnosis.

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© 2020 PriMed Instruments Inc.
Enzol® est une marque déposée de Johnson & Johnson Inc.